LADING- en PAKKETCONTROLE in STOOM STERILISATIE Routine monitoring

4 juni 2015

Ingenious Solutions

Transforming Health

Dirk Eraly

Professional Services Specialist & Scientific support IPD

3M Infection Prevention Solutions

3M Confidential © 3M 2013. All Rights Reserved.

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Lead
the Way

3M Infection Prevention Solutions

3M Confidential © 3M 2013. All Rights Reserved.

3M

The preventable proportion of NI: an overview of 30 published reports

- Great potential exists to decrease NI rates, from a minimum reduction of 10% to a maximum effect of 70% (depending on setting, study design, baseline infection rates and type of infection)
- The most important reduction effect was identified for catheter-related bacteraemia
- **Based on these estimates, at least 20% of all NI are (probably) preventable**

Harbarth S, Sax H, Gastmeier P. JHI 2003 54:258-66

Source: Dettenkofer ESCMID SHEA 2009

Lead
the Way

3M Infection Prevention Solutions

3M Confidential © 3M 2013. All Rights Reserved.

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



Nationale prevalentie (n=54): 7.1%

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Scientific Institute of Public Health 2012;
Source: Dr. M. Goossens (NSIH-PPS)

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
According to the standards: NBN EN ISO 17665

- **NBN EN ISO 17665-1:2006(E)**

- Clause 10 Routine monitoring and control:**

- 10.1 Routine monitoring and control shall be performed on each operating cycle.*

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Lead
the Way

3M Infection Prevention Solutions

3M Confidential © 3M 2013. All Rights Reserved.

3M

NBN EN ISO 17665-1 (NPR-CEN-ISO/TS 17665-2:2009 p 6.1.2.)

Pressure and temperature are not sufficient for monitoring

- **As a consequence the differences between the temperature measured at the temperature measurement point and the temperature calculated from the sterilizer chamber pressure using steam table values (see Annex C) may not be adequate to detect the small volumes of air which could concentrate in lumens or enclosed spaces and prevent steam penetration.**
- Under such circumstances adequate air removal and steam penetration should be predicted from data obtained from a steam penetration test and/or a process monitoring device.
- **Vertaling: Druk en temperatuur tabel om de aanwezigheid van lucht te voorspellen en stoompenetratie in lumen of afgesloten ruimtes kan mogelijk niet correct kunnen zijn.**

Lead
the Way

3M Infection Prevention Solutions

3M Confidential © 3M 2013. All Rights Reserved.

3M

Table B.1 — Example of a schedule of tests for validation and periodic testing

Test and Monitoring	Installation qualification	Operational qualification	Performance qualification	Routine test of the sterilizer	Periodic test of the sterilizer	Comments
Type tests and safety checks: pressure vessel electrical plumbing environmental	×					
Biological challenge test pack		×			×	× in otherwise empty chamber
Air removal test		×		×	×	× in otherwise empty chamber
Air leak test						Sterilizer manufacture test
Physical monitors		×	×	×	×	
Monitoring, recording, control		×	×	×	×	
Independent sensor/recorder						Optional
Biological indicators		×	×	×	×	
Chemical indicators		×	×	×	×	

× = tests that should be considered.

NEN-EN-ISO 15882

(en)

Sterilization of health care products - Chemical indicators - Guidance for selection, use and interpretation of results (ISO 15882:2008, IDT)

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References for load monitoring

- **AJIC Journal: Evaluation of a new rapid readout biological indicator for use in 132°C and 135°C vacuum-assisted steam sterilization cycles, Philip M. Schneider AJIC Journal**
- **Clinical benefits of faster biological indicators, Rose Seavy Healthcare Purchasing N Article 2014**
- **Monitoring of steam sterilization processes in the dental office, J.P.C.M. van Doornmalen, A.G.M. Rietmeijer, A.J. Feilzer, K. Kopinga, Central Service 6/2013**
- **Detecting Steam sterilization proces failures : Bis tell a story. Larson,D., Healthcare Purchasing N 2013**
- **Routine Monitoring of the Steam Sterilization proces in Health Care Facilities.Desinfection, Sterilization and Antisepsis: Principles,Practices, Current issues, New research and New Technologies – APIC Conference Proceedings.(Chapter 9) Ph.Scheider, 2010**

Standards, Recommendations and Guidelines

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WFHSS presentatie Praag 2014 O.Willième Afiso

Currently, hospitals do **MOST** of the right things, on **MOST** patients, **MOST** of the time.

The Checklist helps us do **ALL** the right things, on **ALL** patients, **ALL** the time⁽³⁾

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3M

“Safe Surgery Saves Lives”

WHO Challenge impacts on daily CSSD activities

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Objectives

2

- ▶ To develop your understanding of the aims and global elements of the WHO Surgical Checklist
- ▶ To highlight items in connection with sterilization core business
- ▶ To list impact of this process on daily activities between OR & CSSD
- ▶ To confirm the role of CSSD professional in regard to the “Patient safety checklist” concept (applies to any patient for any type of surgery under any circumstance)

Conclusion

27

- At least two checks of the “patient safety checklist” affect directly or indirectly CSSD professionals work and responsibilities
- Sterilization professionals have a role to play at all times to guarantee security and quality of care
- Their work can be envisaged as essential and rewarding in collaboration with their colleague of the operating room

WHO Surgical Safety Checklist effectief in reduceren post-operatieve complicaties

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






Lead
the Way

© 3M 2008. All Rights Reserved.

3M Infection Prevention
© 3M 2013. All Rights Reserved.

Six commercially available Class 6 Chemical Indicators tested against their stated values. By Doornmalen JPCM van, Hermsen RJ, and Kopinga K. Central Service, 6:400-404, 2012

References:

- ❑ **Six commercially available Class 6 Chemical Indicators tested against their stated values. By Doornmalen JPCM van, Hermsen RJ, and Kopinga K. Central Service, 6:400–404, 2012**
- ❑ **Chemical Indicators—It's not Just A Numbers Game, AAMI ST79:2010 Recommendations for Use of Class 6 Emulating Indicators by Steven Kirckof, BS and Martha Young, MS, BS, CSPDT**
- **Performance of Various Steam Sterilization Indicators Under Optimum and Sub-optimum Exposure Conditions by Philip M. Schneider, Robert R. Reich, Steven S. Kirckof, William G. Foltz, American Journal of Infection Control, Vol. 33, No. 5, Sup. 2, June 2005**
- ❑ **Routine Monitoring of the Steam Sterilization Process in Health Care Facilities, Chapter 9 in Disinfection, Sterilization and Antisepsis: Principles, Practices, Current Issues, New Research and New Technologies – APIC 2009 Pre-Conference Proceedings by P.M. Schneider – LexaMed**
- ❑ **This book chapter contains two studies, both sponsored by 3M Health Care, comparing the performance of Class 6 chemical indicators (CIs) with a Class 5 CI and three biological indicators (BIs) in detecting defined failure conditions in steam sterilization cycles.**


Lead
the Way

Multifunctional Emulating Systems
3M Confidential © 3M 2013. All Rights Reserved.

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Unprocessed

Accept

When the color bar moves anywhere into the "ACCEPT" area, all the critical parameters of steam sterilization have been met.

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© 3M 2007 All Rights Reserved.
Copyright © 2013 3M, Inc.